

# Preventative Maintenance Program



We know the success of your research relies on the precise operation of your Percival Scientific chamber(s). That's why now is the perfect time to consider **Percival's Preventative Maintenance Program**. It's designed to keep your chamber(s) running for years with the accurate performance and environmental control you depend on.

**There's no better way to protect your valuable investment.**

Our service team has over 60 years of experience and can detect a potential issue before it causes a major breakdown.

– Ralan Nelson  
Customer Support Manager



# Preventative Maintenance Program

Through this program, our trained professionals provide onsite annual maintenance services to help minimize the need for costly future repairs, including:

## Operational Checks:


- Calibration of relative humidity and temperature
- Fan motors and blowers
- Electrical System

## Visual Inspections and Cleanings:

- Air-cooled condensers
- Water-cooled condensers, if applicable (for an additional fee)
- Humidifier pans and float valves
- Fan blades (as needed)
- Evaporator coils and drain lines (as needed)

## Gauged Tests:

- Refrigeration systems
- Lighting system
- Humidification system (if applicable)
- CO<sub>2</sub> system (if applicable)
- Designed operating range capabilities



 **PERCIVAL**  
Helping You Create Better Science

To learn more about our Preventative Maintenance Program and get an estimate on pricing, email [service@percival-scientific.com](mailto:service@percival-scientific.com) or call **800.695.2743**.